

ASC Timer Quickstart Guide for OlyOrtho Associates

The ASC timer application is a software utility designed in consultation with Olympia Orthopedic partners that can be used to time each step of an operation and view simple analytics derived from timing operations. This quickstart guide covers how to access and utilize the application, as well as how to interpret some of the data the application presents, as well as what may be some frequently-asked questions.

How to access the ASC Timer

The timer application is a web-based application that you interact with via a web browser. It can be accessed by visiting the following link:

<http://asctimer.olyortho.local/timer/>

This link should work when logged into any device -- computer, phone, or tablet -- as long as the device is configured for and logged into the OlyOrtho local network/VPN.

[Home](#) [Surgeons](#) [Operations](#)

Begin a new operation

Operation type:

Date: ☐ [Today](#)

Surgeon:

Detail:

[Begin new operation](#)

Add a new surgeon

First name:

Last name:

Email:

[Add new Surgeon](#)

Add a new operation type

Operation type:

[Add new operation type](#)

V0.1.0-alpha

© 2025 Nexus Apps LLC. Find a bug? Have a feature request? [Let us know!](#)

[Back to homepage](#)

Standard workflow (with example)

The ASC Timer application provides a flexible way to document and time each step of any arbitrary operation you may perform.

1. Add your information via the `Add a new surgeon` column, and then the `Add a new operation type` to describe the operation type you'd like to time (e.g. 'Total Knee').
2. Fill out the `Begin a new operation` column to kickstart the creation of a new operation instance. You can click/tap the `Today` button to quickly fill out the current date. You can use the detail box to put any information you'd like (such as case numbers) that you'd like associated with this operation.
3. Upon clicking the `Begin new operation` button, you will be redirected to a page where you can submit the steps for the operation you're about to conduct. So for `Total Knee` your step list might look something like this:

Home Surgeons Operations

Edit steps for new operation Total Knee

block done	Remove
tubed	Remove
prepped/draped	Remove
time-out	Remove
incision	Remove
parts in	Remove
closure end	Remove
out of room	Remove
Add Another Step	Submit

V0.1.0-alpha

© 2025 Nexus Apps LLC. Find a bug? Have a feature request? [Let us know!](#)

[Back to homepage](#)

4. upon clicking Submit, you'll be redirected to the timing page. This page is where the timing of the operation occurs.

[Home](#) [Surgeons](#) [Operations](#)

Mark times for Total Knee

In-room

block done: [Now](#)

tubed: [Now](#)

prepped/draped: [Now](#)

time-out: [Now](#)

incision: [Now](#)

parts in: [Now](#)

closure end: [Now](#)

out of room: [Now](#)

Submit

V0.1.0-alpha

© 2025 Nexus Apps LLC. Find a bug? Have a feature request? [Let us know!](#)

[Back to homepage](#)

To start operation timing, click/tap the green `In-room` button when appropriate. Then, as each step of the operation is performed, enter the time that the step is completed (or simply

click/tap the **Now** button to get the current time and fill it in).

[Home](#) [Surgeons](#) [Operations](#)

Mark times for Total Knee

Operation started!

block done: 09:07:21 AM  [Now](#)

tubed: 10:07:24 AM  [Now](#)

prepped/draped: 10:30:30 AM  [Now](#)

time-out: 11:00:37 AM  [Now](#)

incision: 11:23:47 AM  [Now](#)

parts in: 12:07:57 PM  [Now](#)

closure end: 12:45:02 PM  [Now](#)

out of room: 01:13:00 PM  [Now](#)

V0.1.0-alpha

© 2025 Nexus Apps LLC. Find a bug? Have a feature request? [Let us know!](#)

[Back to homepage](#)

5. Upon clicking submit (effectively concluding the timing session), you'll be redirected to the post-operation statistics page.

Operation Instance Total Knee

Case number [here] ... other details

Operation type	Date	Surgeon	Elapsed time
Total Knee	04/22/2025	Tom Foolery	04:02:47

Operation details

Step	Time started	Time ended	Elapsed time	Dist from average
block done	09:00 AM	09:07 AM	00:07:21	79.53% lower than average
tubed	09:07 AM	10:07 AM	01:00:03	28.59% higher than average
prepped/draped	10:07 AM	10:30 AM	00:23:06	53.90% lower than average
time-out	10:30 AM	11:00 AM	00:30:07	34.92% lower than average
incision	11:00 AM	11:23 AM	00:23:10	41.93% lower than average
parts in	11:23 AM	12:07 PM	00:44:10	16.36% lower than average
closure end	12:07 PM	12:45 PM	00:37:05	42.76% lower than average
out of room	12:45 PM	01:13 PM	00:27:58	40.49% lower than average

[Download CSV of this operation](#)

[Update operation instance](#)

[Delete operation instance](#)

V0.1.0-alpha

© 2025 Nexus Apps LLC. Find a bug? Have a feature request? [Let us know!](#)

[Back to homepage](#)

Interpreting the statistics

The `Dist from average` column describes how much more quickly or slowly the corresponding step took compared to other steps of the corresponding type. This calculation will only consider elapsed times of steps performed by the same surgeon, for the same operation type. In other words, the elapsed step times of other surgeons performing the same operation will not factor into the calculation for your step averages. However, if steps vary slightly between operations of the same operation type, or if the order of the steps changes, the elapsed time of those steps is still factored into the calculation.

Frequently-asked questions

***Question:** My browser says 'not secure' when I connect to this application. Is this correct?*

Answer: Currently, the application does not provide the encryption mechanisms required to support HTTPS and therefore display the 'secure' message. This is because the app is only

visible while logged into the OlyOrtho network. HTTPS can be supported in a future update if requested.

***Question:** Can others view my data?*

Answer: Currently, all users of this application can see times and data of all other users. The notion of accounts can be included/supported in a future update if requested.

***Question:** can I access this application from any computer?*

Answer: Yes, the application and any data you save can be accessed from any other computer, phone or tablet that has access to the OlyOrtho local network.

***Question:** I was working through the operation but my tablet died (and/or I accidentally clicked away from the input page). How do I go back?*

Answer: You can view incomplete operation entries (if they exist) by clicking "surgeons" at the top of the application, and then selecting your surgeon entry. You'll find any incomplete entries listed there and will be redirected to the input step you were on before accidentally navigating away from the page.

Home Surgeons Operations			
Tom Foolery			
Complete entries			
Operation type	Date	Surgeon	Elapsed time
Total Knee	03/13/2025	Tom Foolery	07:05:23
Sixth Finger Augmentation	03/13/2025	Tom Foolery	06:04:03
Total Knee	03/13/2025	Tom Foolery	07:19:10
Sixth Finger Augmentation	03/13/2025	Tom Foolery	06:50:35
Total Knee	03/13/2025	Tom Foolery	07:14:21
Sixth Finger Augmentation	03/13/2025	Tom Foolery	05:44:01
Sixth Finger Augmentation	03/13/2025	Tom Foolery	06:14:59
Total Knee	03/13/2025	Tom Foolery	06:05:29
Total Knee	04/22/2025	Tom Foolery	04:02:47
Incomplete entries			
Operation type	Date	Surgeon	Elapsed time
Total Knee	04/22/2025	Tom Foolery	—
Update surgeon entry Delete surgeon entry			

***Question:** what is Nexus Apps LLC? I've never heard of that.*

Answer: Formally, Nexus Apps LLC is a single-member LLC that engages in the engineering and developing of tailored software application solutions. Informally, it's the side gig of a brother-in-law of one of OlyOrtho's partners. :)

***Question:** I've got feedback for the app, and/or an idea about how to do something better, and/or a feature request, and/or I discovered a bug. How can I get this addressed?*

Answer: Any and all feedback is welcome! please email nexusappswa@gmail.com with your feedback and/or ideas and/or feature requests and/or bugs.